


Anti-Pdcd4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q53EI6

Item Number: CM0021069



Introduction

Rabbit polyclonal antibody against human PDCD4, validated for Western blot, IF/ICC, and ELISA. Detects human, mouse, and rat PDCD4. Useful for apoptosis, translation regulation, and tumor suppressor studies.

[Learn More](#)

Parameter	Value
Product Name	PDCD4 Rabbit pAb
Remarks/Alias	H731; PDCD4
Species	Human
Gene ID (Human)	27250
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PDCD4 (NP_055271.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q53EI6
Protein Weight	52kDa
Shipping	Ice bag