


Anti-Pdcd4 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Ihc-P, Elisa - Q53EI6

Item Number: CM0025121



Introduction

KO validated PDCD4 rabbit polyclonal antibody detects human/mouse PDCD4 (SWISS Q53EL6, 52kDa) in WB, IHC-P, ELISA. Targets tumor suppressor inhibiting translation and neoplastic transformation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PDCD4 Rabbit pAb
Remarks/Alias	H731; D4
Species	Human
Gene ID (Human)	27250
Gene ID	27250
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PDCD4 (NP_055271.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q53EL6
Protein Weight	52kDa
Shipping	Ice bag