


Anti-Pdcd5 Polyclonal Antibody For Wb, If/Icc, Elisa - O14737

Item Number: CM0024919



Introduction

High-quality rabbit polyclonal antibody against human PDCD5 (TFAR19). Validated for WB, IF/ICC, ELISA with broad species reactivity (human, mouse, rat). Supports apoptosis and DNA damage response studies.

[Learn More](#)

Parameter	Value
Product Name	PDCD5 Rabbit pAb
Remarks/Alias	TFAR19; PDCD5
Species	Human
GeneID (human)	9141
GeneID	9141
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-125 of human PDCD5 (NP_004699.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14737
Protein Weight	14kDa
Shipping	Ice bag