

Anti-Pdcd7 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q8N8D1

Item Number: CM0003913



Introduction

Recombinant rabbit monoclonal antibody targeting human PDCD7. Validated for western blot, IHC-P, and ELISA; detects PDCD7 in human, mouse, and rat samples. Use in apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	PDCD7 Rabbit mAb
Remarks/Alias	59K, ES18, HES18, PDCD7
Species	Human
Gene ID (Human)	10081
Gene ID	10081
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PDCD7 (Q8N8D1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N8D1
Protein Weight	55kDa
Shipping	Ice bag