


Anti-Pde10A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y233

Item Number: CM0022430



Introduction

Rabbit polyclonal antibody against human PDE10A, a striatal cAMP/cGMP phosphodiesterase involved in signal transduction. Validated for WB, IF/ICC, ELISA; cross-reacts with human and rat. Suitable for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	PDE10A Rabbit pAb
Remarks/Alias	ADSD2, IOLOD, HSPDE10A, PDE10A19, LINC00473, PDE10A
Species	Human
Gene ID (Human)	10846
Gene ID	10846
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-200 of human PDE10A (NP_001124162.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9Y233
Protein Weight	88kDa
Shipping	Ice bag