


Anti-Pde12 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q6L8Q7

Item Number: CM0020454



Introduction

Rabbit polyclonal antibody targeting human PDE12 (2',5'-phosphodiesterase 12) for WB, IF/ICC, ELISA. Reacts with mouse and rat. Suitable for studying antiviral 2-5A system regulation.

[Learn More](#)

Parameter	Specification
Product Name	PDE12 Rabbit pAb
Remarks/Alias	2-PDE; 2'-PDE; PDE12
Species	Human
Gene ID (Human)	201626
Gene ID	201626
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-300 of human PDE12 (NP_001309106.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q6L8Q7
Protein Weight	67kDa
Shipping	Ice bag