

Anti-Pde1A Polyclonal Antibody For Wb And Elisa - P54750

Item Number: CM0013188

Introduction

PDE1A Rabbit pAb is a polyclonal antibody validated for WB and ELISA in human, mouse, and rat. Targets calcium/calmodulin-dependent phosphodiesterase 1A, a key regulator of cyclic nucleotide signaling. Ideal for cardiovascular and neurological research.

[Learn More](#)

Product Name	PDE1A Rabbit pAb
Remarks/Alias	HCAM1; HCAM-1; HSPDE1A; CAM-PDE 1A; CAM-PDE-1A; PDE1A
Species	Human
GeneID (Human)	5136
GeneID	5136
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 446-545 of human PDE1A (NP_005010.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P54750
Protein Weight	61kDa
Shipping	Ice bag