


Anti-Pde1B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q01064

Item Number: CM0022371



Introduction

High-quality rabbit polyclonal antibody targeting human PDE1B, validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat PDE1B, making it ideal for cAMP/cGMP signaling research.

[Learn More](#)

Parameter	Value
Product Name	PDE1B Rabbit pAb
Remarks/Alias	PDE1B1; PDES1B; HEL-S-79p; PDE1B
Species	Human
GeneID (Human)	5153
GeneID	5153
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-536 of human PDE1B (NP_000915.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01064
Protein Weight	61kDa
Shipping	Ice bag