


Anti-Pde3A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14432

Item Number: CM0023920



Introduction

High-quality rabbit polyclonal antibody targeting PDE3A, validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Ideal for apoptosis and cAMP/cGMP signaling research.

[Learn More](#)

Parameter	Value
Product Name	PDE3A Rabbit pAb
Remarks/Alias	HTNB; CGI-PDE; CGI-PDE A; CGI-PDE-A; PDE3A
Species	Human
GeneID (Human)	5139
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 380-660 of human PDE3A (NP_000912.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14432
Protein Weight	125kDa
Shipping	Ice bag