

Anti-Pde3A Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9Z0X4

Item Number: CM0029654

Introduction

Rabbit polyclonal antibody against mouse PDE3A, raised against a recombinant fragment (aa 550-700). Validated for WB, IP, ELISA. Specific for mouse PDE3A, a cGMP-inhibited cAMP phosphodiesterase involved in cyclic nucleotide signaling and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	PDE3A Rabbit pAb
Remarks/Alias	HTNB; CGI-PDE; CGI-PDE A; CGI-PDE-A
Species	Mouse
GeneID (Human)	5139
GeneID	54611
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-700 of mouse PDE3A (NP_061249.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Mouse
SWISS	Q9Z0X4
Protein Weight	125 kDa
Shipping	Ice bag