


Anti-Pde4A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P27815

Item Number: CM0019432



Introduction

Rabbit polyclonal antibody to PDE4A (3',5'-cyclic-AMP phosphodiesterase 4A), validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for cAMP signaling studies and drug screening.

[Learn More](#)

Parameter	Value
Product Name	PDE4A Rabbit pAb
Remarks/Alias	PDE4; DPDE2; PDE46; PDE4A
Species	Human
Gene ID (Human)	5141
Gene ID	5141
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 325-410 of human PDE4A (NP_001104777.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27815
Protein Weight	98kDa
Shipping	Ice bag