

Anti-Pde4B Polyclonal Antibody For Wb, Elisa - Q07343

Item Number: CM0021243



Introduction

High-quality rabbit polyclonal antibody against human PDE4B (DPDE4, PDEIVB). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for cAMP signaling research and drug development studies.

[Learn More](#)

Parameter	Value
Product Name	PDE4B Rabbit pAb
Remarks/Alias	DPDE4, PDEIVB, PDE4B
Species	Human
Gene ID (Human)	5142
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 636-736 of human PDE4B (NP_002591.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q07343
Protein Weight	83kDa
Shipping	Ice bag