

Anti-Pde4C Polyclonal Antibody For If/Icc, Elisa - Q08493

Item Number: CM0017898



Introduction

Rabbit polyclonal antibody targeting human PDE4C (cAMP-specific phosphodiesterase 4C). Validated for IF/ICC and ELISA. Suitable for cAMP signaling studies in human and rat samples.

[Learn More](#)

Parameter	Details
Product Name	PDE4C Rabbit pAb
Remarks/Alias	DPDE1, PDE21, PDE4C
Species	Human
Gene ID (Human)	5143
Gene ID	5143
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 638-712 of human PDE4C (NP_000914.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q08493
Protein Weight	80kDa
Shipping	Ice bag