

Anti-Pde4D Polyclonal Antibody For Wb, Elisa - Q08499

Item Number: CM0018722



Introduction

Rabbit polyclonal antibody targeting human PDE4D, a cAMP-specific phosphodiesterase. Validated for WB and ELISA; cross-reacts with human and rat. Ideal for studying cAMP signaling and related research applications.

[Learn More](#)

Parameter	Value
Product Name	PDE4D Rabbit pAb
Remarks/Alias	DPDE3; PDE43; STRK1; ACRDYS2; HSPDE4D; PDE4DN2; PDE4D
Species	Human
GeneID (Human)	5144
GeneID	5144
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 249-518 of human PDE4D (NP_001184152.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q08499
Protein Weight	91kDa
Shipping	Ice bag