

Anti-Pde4Dip Polyclonal Antibody For Wb, Elisa - Q5Vu43

Item Number: CM0028988

Introduction

Rabbit polyclonal antibody targeting human PDE4DIP (Myomegalin). Validated for WB, ELISA; cross-reacts with mouse. Ideal for microtubule dynamics, cardiac muscle biology, and centrosome research.

[Learn More](#)

Parameter	Value
Product Name	PDE4DIP Rabbit pAb
Remarks/Alias	MMGL; CMYA2; PDE4DIP
Species	Human
Gene ID (Human)	9659
Gene ID	9659
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human PDE4DIP (NP_071754.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q5VU43
Protein Weight	265kDa
Shipping	Ice bag