


Anti-Pde5A Rabbit Polyclonal Antibody For Wb, Elisa - 076074

Item Number: CM0028106



Introduction

Rabbit polyclonal antibody against PDE5A (cGMP-specific phosphodiesterase). Validated for WB and ELISA in human, mouse, and rat. Used in studying cGMP signaling and nitric oxide pathways.

[Learn More](#)

Parameter	Value
Product Name	PDE5A Rabbit pAb
Remarks/Alias	CN5A; PDE5; CGB-PDE; PDE5A
Species	Human
GeneID (Human)	8654
GeneID	8654
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human PDE5A (NP_001074.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	076074
Protein Weight	100kDa
Shipping	Ice bag