

Anti-Pde6A Rabbit Polyclonal Antibody For Wb, Elisa - P16499

Item Number: CM0027185

Introduction

Rabbit polyclonal PDE6A antibody (PDEA, RP43) for WB and ELISA, reactive with human, mouse, rat. Targets rod cGMP phosphodiesterase alpha subunit in visual transduction and retinitis pigmentosa research.

[Learn More](#)

Parameter	Value
Product Name	PDE6A Rabbit pAb
Remarks/Alias	PDEA; RP43; CGPR-A; PDE6A
Species	Human
Gene ID (Human)	5145
Gene ID	5145
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16499
Protein Weight	100kDa
Shipping	Ice bag