

Anti-Pde6B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P35913

Item Number: CM0027841

Introduction

High-quality anti-PDE6B rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, rat PDE6B, a key enzyme in phototransduction. Ideal for vision research.

[Learn More](#)

Parameter	Value
Product Name	PDE6B Rabbit pAb
Remarks/Alias	rd1; PDEB; RP40; CSNB3; CSNBAD2; GMP-PDEbeta; PDE6B
Species	Human
GeneID (Human)	5158
GeneID	5158
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human PDE6B (NP_001138763.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35913
Protein Weight	98kDa
Shipping	Ice bag