

Anti-Pde6D Polyclonal Antibody For Wb, If/Icc, Elisa - O43924

Item Number: CM0016962



Introduction

CamelBio's PDE6D Polyclonal Antibody (Rabbit pAb) validated for WB, IF/ICC, ELISA in human, mouse, rat. Ideal for studying prenylated protein trafficking and ciliary targeting in research and diagnostics. One-stop IVD raw material sourcing.

[Learn More](#)

Parameter	Value
Product Name	PDE6D Rabbit pAb
Remarks/Alias	PDED; JBTS22; PDE6D
Species	Human
Gene ID (Human)	5147
Gene ID	5147
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human PDE6D (NP_002592.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43924
Protein Weight	17kDa
Shipping	Ice bag