

Anti-Pde7A Polyclonal Antibody For Wb, Elisa - Q13946

Item Number: CM0016459



Introduction

High-affinity cAMP-specific phosphodiesterase 7A rabbit polyclonal antibody validated for WB and ELISA applications in human, mouse, and rat samples. Ideal for studying cAMP signaling and muscle biology.

[Learn More](#)

Parameter	Value
Product Name	PDE7A Rabbit pAb
Remarks/Alias	HCP1; PDE7; PDE7A
Species	Human
Gene ID (Human)	5150
Gene ID	5150
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 345-456 of human PDE7A (NP_002594.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13946
Protein Weight	56kDa
Shipping	Ice bag