


Anti-Pde8A Rabbit Polyclonal Antibody For Wb, Elisa - O60658

Item Number: CM0017232



Introduction

Rabbit polyclonal antibody targeting human PDE8A, validated for WB and ELISA with mouse and rat cross-reactivity. Suitable for cAMP signaling and cell death research.

[Learn More](#)

Parameter	Value
Product Name	PDE8A Rabbit pAb
Remarks/Alias	HsT19550, PDE8A
Species	Human
GeneID (human)	5151
GeneID	5151
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human PDE8A (NP_002596.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O60658
Protein Weight	93kDa
Shipping	Ice bag