


Anti-Pdf Rabbit Polyclonal Antibody For Wb And Elisa - Q9Hbh1

Item Number: CM0018048



Introduction

Anti-PDF rabbit polyclonal antibody validated for WB and ELISA in human, mouse, rat. Targets mitochondrial peptide deformylase for protein biosynthesis studies and IVD development.

[Learn More](#)

Parameter	Value
Product Name	PDF Rabbit pAb
Remarks/Alias	PDF
Species	Human
GeneID	64146
Immunogen	Recombinant fusion protein containing amino acids 40-243 of human PDF (NP_071736.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HBH1
Protein Weight	27 kDa
Shipping	Ice bag