


Anti-Pdgfa Polyclonal Antibody For Ihc-P, Elisa - P04085

Item Number: CM0021387



Introduction

Rabbit polyclonal antibody against human PDGFA, validated for IHC-P and ELISA. Cross-reacts with mouse. Suitable for developmental, cell proliferation, and wound healing research.

[Learn More](#)

Parameter	Value
Product Name	PDGFA Rabbit pAb
Remarks/Alias	PDGF1; PDGF-A; PDGFA
Species	Human
GeneID (Human)	5154
GeneID	5154
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-110 of human PDGFA (NP_002598.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P04085
Protein Weight	24kDa
Shipping	Ice bag