

Anti-Pdgfb Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P01127

Item Number: CM0016764



Introduction

Polyclonal antibody targeting human PDGFB (platelet-derived growth factor subunit B), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with rat; detects a 27 kDa band. Useful for cell proliferation, angiogenesis, and developmental biology studies.

[Learn More](#)

Parameter	Specification
Product Name	PDGFB Rabbit pAb
Remarks/Alias	SIS; SSV; IBGC5; PDGF2; c-sis; PDGF-2; PDGFB
Species	Human
Gene ID (Human)	5155
Gene ID	5155
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 140-241 of human PDGFB (NP_002599.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P01127
Protein Weight	27kDa
Shipping	Ice bag