


Anti-Pdgfb Rabbit Monoclonal Antibody For Wb, Ihc-P, If-P, Elisa - P01127

Item Number: CM0007606



Introduction

Rabbit monoclonal antibody against human PDGFB (P01127). Suitable for WB, IHC-P, IF-P, ELISA. Cross-reacts with mouse and rat. Validated with synthetic peptide immunogen. Ideal for cell proliferation and oncology studies.

[Learn More](#)

Parameter	Value
Product Name	PDGFB Rabbit mAb
Remarks/Alias	SIS; SSV; IBGC5; PDGF2; c-sis; PDGF-2; PDGFB
Species	Human
GeneID	5155
Immunogen	Synthetic peptide corresponding to amino acids 142-241 of human PDGFB (NP_002599.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01127
Protein Weight	27 kDa
Shipping	Ice bag