

Anti-Pdgfc Polyclonal Antibody For Wb, Elisa - Q9Nra1

Item Number: CM0018205

Introduction

Rabbit polyclonal antibody against human PDGFC, validated for WB and ELISA. Detects human, mouse, and rat PDGFC, a growth factor involved in embryonic development, angiogenesis, and wound healing.

[Learn More](#)

Parameter	Value
Product Name	PDGFC Rabbit pAb
Remarks/Alias	SCDGF; FALLOTEIN; PDGFC
Species	Human
Gene ID (Human)	56034
Gene ID	56034
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 166-345 of human PDGFC (NP_057289.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRA1
Protein Weight	39kDa
Shipping	Ice bag