

Anti-Pdgfd Rabbit Polyclonal Antibody For Wb And Elisa - Q9Gzp0

Item Number: CM0019720

Introduction

Rabbit polyclonal antibody targeting human PDGFD (IEGF, SCDGF-B). Validated for Western blot and ELISA, with cross-reactivity to mouse. Suitable for research on cell proliferation, angiogenesis, wound healing, and oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	PDGFD Rabbit pAb
Remarks/Alias	IEGF; SCDGFB; MSTP036; SCDGF-B; PDGFD
Species	Human
Gene ID (Human)	80310
Gene ID	80310
Immunogen	Recombinant protein corresponding to a sequence within amino acids 19-307 of human PDGFD (NP_079484.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9GZP0
Protein Weight	43kDa
Shipping	Ice bag