

Anti-Pdgfd Rabbit Polyclonal Antibody For Wb, Elisa - Q9Gzp0

Item Number: CM0025084

Introduction

Anti-PDGFD rabbit polyclonal antibody validated for WB and ELISA with human and mouse reactivity. Targets the platelet-derived growth factor D involved in cell proliferation, migration, and angiogenesis. Suitable for cancer and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	PDGFD Rabbit pAb
Remarks/Alias	IEGF; SCDGFB; MSTP036; SCDGF-B; PDGFD
Species	Human
Gene ID	80310
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human PDGFD (NP_079484.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9GZP0
Protein Weight	43kDa
Shipping	Ice bag