

Anti-Pdgfd Rabbit Polyclonal Antibody For Wb, Elisa - Q9Gzp0

Item Number: CM0013058

Introduction

PDGFD Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse, and rat PDGFD. Ideal for research on cell proliferation, development, and wound healing. Q9GZP0.

[Learn More](#)

Parameter	Value
Product Name	PDGFD Rabbit pAb
Remarks/Alias	IEGF; SCDGFB; MSTP036; SCDGF-B; PDGFD
Species	Human
GeneID (human)	80310
GeneID	80310
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZP0
Protein Weight	43kDa
Shipping	Ice bag