

Anti-Pdgfra/Cd140A Rabbit Monoclonal Antibody () For Flow Cytometry - P16234

Item Number: CM0009595



Introduction

Recombinant rabbit monoclonal antibody conjugated to , targeting human and mouse PDGFR α /CD140a. Validated for flow cytometry, ideal for receptor tyrosine kinase signaling and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Mouse PDGFR α /CD140a mAb
Remarks/Alias	CD140A; PDGFR2; PDGFR-2
Species	Human
GeneID	5156
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-524 of human PDGFR α (NP_006197.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	P16234
Protein Weight	24kDa/82kDa/123kDa
Shipping	Ice bag