

Anti-Pdgfra/Pdgfr Alpha Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P16234

Item Number: CM0019715



Introduction

High-quality rabbit polyclonal antibody against human PDGFR alpha (CD140a). Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for research on receptor tyrosine kinase signaling, cancer, and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	PDGFR alpha Rabbit pAb
Remarks/Alias	CD140A; PDGFR2; PDGFR-2; PDGFR alpha
Species	Human
GeneID (human)	5156
GeneID	5156
Immunogen	Recombinant protein, recombinant fusion protein containing a sequence corresponding to amino acids 896-1059 of human PDGFR alpha (NP_006197.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16234
Protein Weight	123kDa
Shipping	Ice bag