


Anti-Pdgfrb Polyclonal Antibody For Wb - P09619

Item Number: CM0017703



Introduction

Rabbit polyclonal antibody against PDGFRB (P09619) for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for IVD research and development.

[Learn More](#)

Parameter	Value
Product Name	PDGFRB Rabbit pAb
Remarks/Alias	IMF1; KOGS; IBGC4; JTK12; PDGFR; PENTT; CD140B; PDGFR1; PDGFR-1; PDGFRB
Species	Human
GeneID (Human)	5159
Gene ID	5159
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 900-1000 of human PDGFRB (NP_002600.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09619
Protein Weight	124 kDa
Shipping	Ice bag