


Anti-Pdgfrb Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P09619

Item Number: CM0002852



Introduction

Rabbit monoclonal antibody against human PDGFRB (CD140B). Validated for Western blot, IF/ICC, and ELISA. High specificity for PDGFRB; cross-reacts with human and mouse. Suitable for studying PDGF signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	PDGFRB Rabbit mAb
Remarks/Alias	IMF1; KOGS; IBGC4; JTK12; PDGFR; PENTT; CD140B; PDGFR1; PDGFR-1; PDGFRB
Species	Human
GeneID (Human)	5159
GeneID	5159
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1000-1106 of human PDGFRB (P09619).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09619
Protein Weight	124kDa
Shipping	Ice bag