

Anti-Pdgfrl Rabbit Polyclonal Antibody For Wb, Elisa - Q15198

Item Number: CM0027660

Introduction

PDGFRL (PDGRL/PRLTS) rabbit polyclonal antibody detects human and mouse platelet-derived growth factor receptor-like protein, a tumor suppressor. Suitable for WB and ELISA, 42 kDa.

[Learn More](#)

Parameter	Value
Product Name	PDGFRL Rabbit pAb
Remarks/Alias	PDGRL; PRLTS; PDGFRL
Species	Human
GeneID (Human)	5157
GeneID	5157
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 15-205 of human PDGFRL (NP_006198.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15198
Protein Weight	42kDa
Shipping	Ice bag