

Anti-Pdha1 Monoclonal Antibody For Wb - P08559

Item Number: CM0004488



Introduction

Rabbit monoclonal antibody specific to human PDHA1, suitable for WB and ELISA. Cross-reacts with mouse PDHA1. Targets a key mitochondrial enzyme linking glycolysis and TCA cycle.

[Learn More](#)

Parameter	Value
Product Name	PDHA1 Rabbit mAb
Remarks/Alias	PDHA; PDHAD; PHE1A; E1alpha; PDHCE1A; PDHA1
Species	Human
Gene ID (Human)	5160
Gene ID	5160
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 248-348 of human PDHA1 (NP_000275.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08559
Protein Weight	43kDa
Shipping	Ice bag