

Anti-Pdha1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P08559

Item Number: CM0018799



Introduction

High-quality rabbit polyclonal antibody targeting human PDHA1, validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Ideal for studying mitochondrial pyruvate dehydrogenase complex and metabolism.

[Learn More](#)

Parameter	Value
Product Name	PDHA1 Rabbit pAb
Remarks/Alias	PDHA; PDHAD; PHE1A; E1alpha; PDHCE1A; PDHA1
Species	Human
GeneID (human)	5160
GeneID	5160
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-390 of human PDHA1 (NP_000275.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08559
Protein Weight	43kDa
Shipping	Ice bag