

Anti-Pdha2 Rabbit Polyclonal Antibody For Wb, If/Icc, And Elisa - P35487

Item Number: CM0016013

Introduction

Rabbit polyclonal antibody against PDHA2, validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Detects a 43 kDa mitochondrial protein in testis, suitable for pyruvate dehydrogenase complex studies.

[Learn More](#)

Parameter	Value
Product Name	PDHA2 Rabbit pAb
Remarks/Alias	Pdhal; PDHA2
Species	Mouse
Gene ID (Human)	5161
Gene ID (Mouse)	18598
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 230-391 of mouse PDHA2 (NP_032837.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35487
Protein Weight	43kDa
Shipping	Ice bag