

Anti-Pdhb Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P11177

Item Number: CM0009566



Introduction

Rabbit monoclonal antibody against human PDHB, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for metabolic research and assay development.

[Learn More](#)

Parameter	Value
Product Name	PDHB Rabbit mAb
Remarks/Alias	PDHBD; PHE1B; E1beta; PDHE1B; PDHE1-B; PDHB
Species	Human
Gene ID (Human)	5162
Gene ID	5162
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 260-359 of human Pyruvate Dehydrogenase E1 beta subunit (P11177).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11177
Protein Weight	39kDa
Shipping	Ice bag