

Anti-Pdhb Rabbit Polyclonal Antibody For Wb And Elisa - P11177

Item Number: CM0026079



Introduction

PDHB Rabbit Polyclonal Antibody validated for Western blot and ELISA. Detects human, mouse, rat PDHB. Suitable for research on pyruvate dehydrogenase complex, glucose metabolism, and mitochondrial disorders.

[Learn More](#)

Parameter	Value
Product Name	PDHB Rabbit pAb
Remarks/Alias	PDHBD; PHE1B; E1beta; PDHE1B; PDHE1-B; PDHB
Species	Human
Gene ID (Human)	5162
Gene ID	5162
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PDHB (NP_000916.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11177
Protein Weight	39kDa
Shipping	Ice bag