

Anti-Pdia2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q13087

Item Number: CM0015289

Introduction

Rabbit polyclonal antibody against human PDIA2, suitable for WB, IF/ICC, and ELISA. Detects human, mouse, and rat targets. Useful for studying endoplasmic reticulum chaperone and estrogen-binding functions.

[Learn More](#)

Parameter	Value
Product Name	PDIA2 Rabbit pAb
Remarks/Alias	PDI; PDA2; PDIP; PDIR; PDIA2
Species	Human
GeneID (Human)	64714
GeneID	64714
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 326-525 of human PDIA2 (NP_006840.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13087
Protein Weight	58kDa
Shipping	Ice bag