

Anti-Pdia4 Rabbit Polyclonal Antibody - P13667

Item Number: CM0027516



Introduction

Rabbit polyclonal antibody targeting human PDIA4 (ERP70/ERP72). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for protein disulfide isomerase and ER stress research.

[Learn More](#)

Parameter	Value
Product Name	PDIA4 Rabbit pAb
Remarks/Alias	ERP70; ERP72; ERp-72; PDIA4
Species	Human
GeneID (Human)	9601
GeneID	9601
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human PDIA4 (NP_004902.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13667
Protein Weight	73kDa
Shipping	Ice bag