

Anti-Pdia6 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q15084

Item Number: CM0009676



Introduction

High-quality rabbit monoclonal antibody targeting human PDIA6 (ERP5/P5), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat; suitable for studying ER stress, unfolded protein response, and platelet function.

[Learn More](#)

Parameter	Value
Product Name	PDIA6 Rabbit mAb
Remarks/Alias	P5, ERP5, TXNDC7, PDIA6
Species	Human
Gene ID	10130
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PDIA6 (NP_005733.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15084
Protein Weight	48kDa
Shipping	Ice bag