

Anti-Pdia6 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15084

Item Number: CM0024829



Introduction

High-quality rabbit polyclonal antibody targeting human PDIA6 (ERP5/P5). Validated for WB, IHC-P, ELISA; cross-reacts with mouse & rat. Ideal for platelet and ER stress research. ~48 kDa.

[Learn More](#)

Parameter	Value
Product Name	PDIA6 Rabbit pAb
Remarks/Alias	P5, ERP5, TXNDC7, PDIA6
Species	Human
GeneID (Human)	10130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human PDIA6 (NP_005733.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15084
Protein Weight	48kDa
Shipping	Ice bag