

Anti-Pdk1/Pdhk1 Rabbit Monoclonal Antibody For Wb, Elisa - Q15118

Item Number: CM0004948



Introduction

Recombinant rabbit monoclonal antibody against human PDK1/PDHK1 (Q15118), validated for WB and ELISA. Cross-reacts with mouse. Ideal for studying glucose metabolism and hypoxia signaling.

[Learn More](#)

Parameter	Value
Product Name	PDK1/PDHK1 Rabbit mAb
Remarks/Alias	PDK1; pyruvate dehydrogenase kinase 1; PDK1/PDHK1
Species	Human
GeneID (Human)	5163
GeneID	5163
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PDK1/PDHK1 (Q15118).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15118
Protein Weight	49kDa
Shipping	Ice bag