

Anti-Pdk1/Pdhk1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q15118

Item Number: CM0008755



Introduction

Rabbit monoclonal antibody against PDK1/PDHK1, validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studying glucose metabolism, hypoxia, and mitochondrial function.

[Learn More](#)

Parameter	Value
Product Name	PDK1/PDHK1 Rabbit mAb
Remarks/Alias	PDK1; pyruvate dehydrogenase kinase 1; PDK1/PDHK1
Species	Human
Gene ID (Human)	5163
Gene ID	5163
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15118
Protein Weight	49kDa
Shipping	Ice bag