

Anti-Pdk1/Pdhk1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15118

Item Number: CM0015701



Introduction

High-quality rabbit polyclonal antibody against human PDK1/PDHK1 (pyruvate dehydrogenase kinase 1). Validated for WB, IF/ICC, and ELISA, with cross-reactivity to mouse and rat. Ideal for studying glucose metabolism regulation and hypoxia response. ~49 kDa.

[Learn More](#)

Parameter	Value
Product Name	PDK1/PDHK1 Rabbit pAb
Remarks/Alias	PDK1; PDK1/PDHK1
Species	Human
Gene ID (Human)	5163
Gene ID	5163
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 337-436 of human PDK1/PDHK1 (NP_002601.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15118
Protein Weight	49kDa
Shipping	Ice bag