


Anti-Pdk1/Pdpk1 Rabbit Polyclonal Antibody For Wb And Elisa - O15530

Item Number: CM0023024



Introduction

Rabbit polyclonal antibody against PDK1/PDPK1, suitable for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying insulin signaling, cell proliferation, and survival pathways.

[Learn More](#)

Parameter	Value
Product Name	PDK1/PDPK1 Rabbit pAb
Remarks/Alias	PDK1, PDPK2, PDPK2P, PRO0461, PDK1/PDPK1
Species	Human
Gene ID (Human)	5170
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 409-556 of human PDK1/PDPK1 (NP_002604.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15530
Protein Weight	63kDa
Shipping	Ice bag