

Anti-Pdk2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q15119

Item Number: CM0012657



Introduction

High-specificity rabbit monoclonal antibody targeting human PDK2 (pyruvate dehydrogenase kinase 2). Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for metabolism and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	PDK2 Rabbit mAb
Remarks/Alias	PDHK2; PDKII; PDK2
Species	Human
GeneID (Human)	5164
GeneID	5164
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 308-407 of human PDK2 (Q15119).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15119
Protein Weight	39kDa/46kDa
Shipping	Ice bag