

Anti-Pdk2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15119

Item Number: CM0026357



Introduction

Anti-PDK2 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, and rat PDK2. Ideal for studying glucose and fatty acid metabolism regulation.

[Learn More](#)

Parameter	Value
Product Name	PDK2 Rabbit pAb
Remarks/Alias	PDHK2; PDKII; PDK2
Species	Human
GeneID (Human)	5164
GeneID	5164
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 128-407 of human PDK2 (NP_002602.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15119
Protein Weight	46kDa
Shipping	Ice bag