


Anti-Pdk3 Polyclonal Antibody For Wb, Elisa - Q15120

Item Number: CM0030246



Introduction

PDK3 Rabbit pAb detects human, mouse, rat PDK3. Suitable for WB and ELISA, this antibody targets a mitochondrial kinase involved in glucose metabolism regulation.

[Learn More](#)

Parameter	Value
Product Name	PDK3 Rabbit pAb
Remarks/Alias	CMTX6; GS1-358P8.4; PDK3
Species	Human
GeneID (Human)	5165
GeneID	5165
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 137-406 of human PDK3 (NP_005382.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15120
Protein Weight	47kDa
Shipping	Ice bag